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● MAIN ● POLICIES ● SAFETY ● PROCEDURES
● EQUIPMENT ● EXPERIMENTS ● CHARTS/TABLES ● INDEX

1. Safety
2. Objectives/Overview
3. Procedures
4. Observations
5. Equipment
6. View Techniques
7. Waste Disposal/Clean-up
8. Conclusions
9. Calculations/Set-Ups
10. Grading Scale
11. Review Prelab Questions
12. Review Postlab Questions

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OBJECTIVE/SUMMARY

The primary objective of this experiment is to determine the concentration of an unknown nickel (II) sulfate solution. You will be using the Colorimeter shown in Figure 1. In this device, red light from the LED light source will pass through the solution and strike a photocell. The NiSO_4 solution used in this experiment has a deep green color. A higher concentration of the colored solution absorbs more light and transmits less light than a solution of lower concentration. The computer-interfaced colorimeter monitors the light received by the photocell as either an absorbance or a percent transmittance value.

Figure 1

Figure 2

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
● MAIN ● POLICIES ● SAFETY ● PROCEDURES
● EQUIPMENT ● EXPERIMENTS ● CHARTS/TABLES ● INDEX

1. Safety
2. Objectives/Overview
3. Procedures
4. Observations
5. Equipment
6. View Techniques
7. Waste Disposal/Clean-up
8. Conclusions
9. Calculations/Set-Ups
10. Grading Scale
11. Review Prelab Questions
12. Review Postlab Questions

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Prepare Solutions
Keep the remaining 0.400M NiSO_4 in the 100mL beaker to use in the 5th trial.
Approximate volumes and concentrations for the trials are shown left. Actual values of volume should be recorded on the data sheet. These will then be used to calculate actual concentrations.

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Experiment 13

Determining the Concentration of a Solution: Beer's Law

OBJECTIVES/OVERVIEW

The primary objective of this experiment is to determine the concentration of an unknown nickel (II) sulfate solution. You will be using the Colorimeter shown in Figure 1. In this device, red light from the LED light source will pass through the solution and strike a photocell. The NiSO₄ solution used in this experiment has a deep green color. A higher concentration of the colored solution absorbs more light and transmits less light than a solution of lower concentration. The computer-interfaced colorimeter monitors the light received by the photocell as either an absorbance or a percent transmittance value.



Figure 1

Figure 2


8. When you have finished reading and entering all five of your solutions, click on "STOP" under the "EXPERIMENT" menu.

9. Copy the absorbance and concentration data pairs from the Data Table on the screen, into the Data and Calculations Table of your report.

10. Examine the graph of absorbance versus concentration presented on your computer screen. Does it

12. Click "INTERPOLATE" from the "ANALYZE" menu. A vertical cursor will appear on the screen. As you move the cursor back and forth across the screen with your mouse, you will see a small "zero" moving up and down your graph. Look in the information box connected to your graph. You will see changing values. Adjust the position of the cursor until you see the value of the y-axis equal to the value of the absorbance of your unknown. The corresponding value of the x-axis is the concentration of your unknown. Write this value into the appropriate space in your Data and Calculation Table. Leave the "INTERPOLATE" cursor on the screen so it is a part of the graph you print.

13. Select "PRINT" from the "FILE" menu. Enter your name and your partner's name. You may also enter comments. You need to print two copies of the graph, one for each partner. Several parts of your graph will need to be labeled: the calibration line prepared with known concentrations; the interpolation line; and the value of the unknown's concentration.

Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

Experiment 13
Determining the Concentration of a Solution: Beer's Law

1. [Safety](#)
2. [Objectives/Overview](#)
3. [Procedures](#)
4. [Observations](#)
5. [Equipment](#)
6. [View Techniques](#)
7. [Waste Disposal/Clean-up](#)
8. [Conclusions](#)
9. [Calculations/Set-Ups](#)
10. [Grading Scale](#)
11. [Review Prelab Questions](#)
12. [Review Postlab Questions](#)

Prepare Solutions
Dissolve the remaining 0.450M NiSO₄ in the 100mL beaker to use in the titration.
Approximate volumes and concentrations for the trials are shown left. Actual values of volume should be recorded on the data sheet. These will then be used to calculate actual concentrations.
Press forward to continue.

100 mL NiSO ₄	100 mL NiSO ₄	100 mL NiSO ₄	100 mL NiSO ₄	100 mL NiSO ₄
0.450M	0.240M	0.160M	0.080M	0.040M

Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

Experiment 13
Determining the Concentration of a Solution: Beer's Law

1. [Safety](#)
2. [Objectives/Overview](#)
3. [Procedures](#)
4. [Observations](#)
5. [Equipment](#)
6. [View Techniques](#)
7. [Waste Disposal/Clean-up](#)
8. [Conclusions](#)
9. [Calculations/Set-Ups](#)
10. [Grading Scale](#)
11. [Review Prelab Questions](#)
12. [Review Postlab Questions](#)

Calibrate Colorimeter
Always position the cuvette with its smooth sides facing toward the reference mark at the right of the cuvette slot on the colorimeter.
Press forward to continue.

Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

Experiment 13
Determining the Concentration of a Solution: Beer's Law

1. [Safety](#)
2. [Objectives/Overview](#)
3. [Procedures](#)
4. [Observations](#)
5. [Equipment](#)
6. [View Techniques](#)
7. [Waste Disposal/Clean-up](#)
8. [Conclusions](#)
9. [Calculations/Set-Ups](#)
10. [Grading Scale](#)
11. [Review Prelab Questions](#)
12. [Review Postlab Questions](#)

Calibration Curve
You should enter the concentration of your solution in the box (in this case 0.080). Then, press the enter key on the computer keyboard. The data pair you just indicated should now be plotted on the graph and entered in the data table.
Enter '0.080' in the concentration field and click 'Enter' on your keyboard.

Absorbance-Live (A)
0.089

Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

Experiment 13
Determining the Concentration of a Solution: Beer's Law

1. [Safety](#)
2. [Objectives/Overview](#)
3. [Procedures](#)
4. [Observations](#)
5. [Equipment](#)
6. [View Techniques](#)
7. [Waste Disposal/Clean-up](#)
8. [Conclusions](#)
9. [Calculations/Set-Ups](#)
10. [Grading Scale](#)
11. [Review Prelab Questions](#)
12. [Review Postlab Questions](#)

Calibration Curve
Examine the graph of absorbance versus concentration presented on your computer screen. Your graph should look like the one pictured left.
Click forward to continue.

Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

Experiment 13
Determining the Concentration of a Solution: Beer's Law

1. [Safety](#)
2. [Objectives/Overview](#)
3. [Procedures](#)
4. [Observations](#)
5. [Equipment](#)
6. [View Techniques](#)
7. [Waste Disposal/Clean-up](#)
8. [Conclusions](#)
9. [Calculations/Set-Ups](#)
10. [Grading Scale](#)
11. [Review Prelab Questions](#)
12. [Review Postlab Questions](#)

Clean Up
Nickel Sulfate waste should be poured into the labeled waste container. Wash equipment in the sink at your lab station and return the clean, borrowed equipment (cuvettes, beakers, and watch glasses) to their storage; pipettes should be placed in the washing basket, point end up.

Absorbance vs. Concentration

$y = mx + b$
 $m = 0.450$ $b = 0.000$
 $x = 0.078$ $y = 0.0517$

Click [HERE](#) to view animation again.