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Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

● MAIN ● POLICIES ● SAFETY ● PROCEDURES
● EQUIPMENT ● EXPERIMENTS ● CHARTS/TABLES ● INDEX

1. Safety
2. Objectives/Overview
3. Procedures
4. Observations
5. Equipment
6. View Techniques
7. Waste Disposal/Clean-up
8. Conclusions
9. Calculations/Set-Ups
10. Grading Scale
11. Review Prelab Questions
12. Review Postlab Questions

Experiment 13
Determining the Concentration of a Solution: Beer's Law

OBJECTIVE/S/OVERVIEW

The primary objective of this experiment is to determine the concentration of an unknown nickel (II) sulfate solution. You will be using the Colorimeter shown in Figure 1. In this device, red light from the LED light source will pass through the solution and strike a photocell. The NiSO_4 solution used in this experiment has a deep green color. A higher concentration of the colored solution absorbs more light and transmits less light than a solution of lower concentration. The computer-interfaced colorimeter monitors the light received by the photocell as either an absorbance or a percent transmittance value.


Figure 1


Figure 2


Multimedia **Preparation** for Chemistry 1151 and 1161 Laboratory

● MAIN ● POLICIES ● SAFETY ● PROCEDURES
● EQUIPMENT ● EXPERIMENTS ● CHARTS/TABLES ● INDEX

1. Safety
2. Objectives/Overview
3. Procedures
4. Observations
5. Equipment
6. View Techniques
7. Waste Disposal/Clean-up
8. Conclusions
9. Calculations/Set-Ups
10. Grading Scale
11. Review Prelab Questions
12. Review Postlab Questions

Experiment 13
Determining the Concentration of a Solution: Beer's Law

TITLE

15mL NiSO_4 , 6mL H_2O , 0.400M
8mL NiSO_4 , 2mL H_2O , 0.240M
6mL NiSO_4 , 4mL H_2O , 0.160M
4mL NiSO_4 , 6mL H_2O , 0.0800M
2mL NiSO_4 , 8mL H_2O , 0.0800M

Prepare Solutions
Keep the remaining 0.400M NiSO_4 in the 100mL beaker to use in the fifth trial.

Approximate volumes and concentrations for the trials are shown left. Actual values of volume should be recorded on the data sheet. These will then be used to calculate actual concentrations.

Press forward to continue.